



	MED-XN91	MED-XN71	MED-XN51LP
Bus Interface	PCI- Express × 16	PCI-Express × 16	PCI- Express × 16
Compatible OS	Windows 10, Windows 8.1, Windows 7 (four output max)	Windows 10, Windows 8.1, Windows 7 (four output max)	Windows 10, Windows 8.1, Windows 7 (four output max)
Frame Buffer Memory	8 GB	5 GB	4 GB
Display Grayscale Tones/Colors	10 Bit, 8 Bit	10 Bit, 8 Bit	10 Bit, 8 Bit
Output Terminal	4x DisplayPort	4x DisplayPort	4x Mini DisplayPort
Cable	1x cable (DisplayPort – DVI-D)	1x cable (DisplayPort – DVI-D)	2x cable (Mini DisplayPort – DisplayPort) 1x cable (Mini DisplayPort – DVI-D)
Daisy Chain Support	◆	◆	◆
Maximum Power Consumption	105 W	75 W	47 W
Chassis	Standard	Standard	Standard and Low-profile
Dimensions (W × H)	241.3 × 111.2 mm	200.7 × 111.2 mm	153.9 × 68.9 mm
Compatibility			
8MP RX850	★	✓	✓
6MP RX660	✓	★	✓
5MP SVP RX560-MD	★	✓	✓
5MP RX560	★	✓	✓
3MP RX360	✓	★	✓
3MP RX250	✓	✓	★
5MP SVP GX560-MD	★	✓	✓
5MP GX560	★	✓	✓
3MP GX340	✓	★	✓
2MP GX240	✓	✓	★
8MP MX315W	✓	★	✓
23MP MX242W	✓	✓	★
2MP MX216-HB	✓	✓	★
2MP MX216-SB	✓	✓	★
1MP MX194	✓	✓	✓
2MP MS236WT	✓	✓	✓

✓ Compatible ★ Recommended

Austria | www.eizo.at
Belgium | www.eizo.be
Czech Republic | www.eizo.cz
United Kingdom | www.eizo.co.uk

Germany | www.eizo.de
Hungary | www.eizo.hu
Italy | www.eizo.it

The Netherlands | www.eizo.nl
Nordics | www.eizo.se
Slovakia | www.eizomonitor.sk